

TC WIRE TWISTER

straight; serrated

18.5cm (7.4")




SKU: 051-242-18

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Wire twisting Forceps
Aesculap	BM113
MATERIAL	Tungsten Carbide(TC)
KLS Martin	22-506-18
Medicon	14.44.67
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Wire twisting Forceps
Aesculap	BM113
MATERIAL	Tungsten Carbide(TC)

KLS Martin

22-506-18

Medicon

14.44.67

SPECIFICATIONS

[serrated](#), [straight](#)