

TELESCOPE HD+2 FOR LAPAROSCOPY

HYSTEROSCOPY, CYSTOSCOPY

Ø 2.9mm; 30°; 300 mm; autoclavable



SKU: 601-103-30HD2

Categories: [Endoscopy](#), [Equipment](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Telescopes
SCOPE	GUBBINI/Gyn , Gynecology , Laparoscopy , Urology



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Telescopes
SCOPE	GUBBINI/Gyn , Gynecology , Laparoscopy , Urology