

## TENACULUM FORCEPS Ø 1.6MM; FR. 5; 34CM

semi-rigid; oval; double action

Luer-Lock



**SKU:** 583-278-16LL

**Categories:** [Endoscopy](#), [GUBBINI/Gyn](#), [Gynecology](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)  
**RANGE** [Micro Instruments Hysteroscopy](#)  
**TYPE** ["Semi-Rigid" with LuerLock](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)  
**RANGE** [Micro Instruments Hysteroscopy](#)  
**TYPE** ["Semi-Rigid" with LuerLock](#)