

TENACULUM FORCEPS Ø 1.6MM; FR. 5; 34CM

semi-rigid; oval; double action

Luer-Lock



SKU: 583-278-16LL

Categories: [Endoscopy](#), [GUBBINI/Gyn](#), [Gynecology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Micro Instruments Hysteroscopy](#)
TYPE ["Semi-Rigid" with LuerLock](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Micro Instruments Hysteroscopy](#)
TYPE ["Semi-Rigid" with LuerLock](#)