

TENDON HOLDING FORCEPS

straight; with fixation screw

11.5cm (4.6")




SKU: 195-252-11

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Seizing Instruments
Aesculap	BF307
KLS Martin	23-680-13
Medicon	78.12.13
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Seizing Instruments
Aesculap	BF307
KLS Martin	23-680-13
Medicon	78.12.13

SPECIFICATIONS

[straight](#)