

## TENDON HOLDING FORCEPS

straight; with fixation screw

11.5cm (4.6")




**SKU:** 195-252-11

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tendon Seizing Instruments</a>
<b>Aesculap</b>	BF307
<b>KLS Martin</b>	23-680-13
<b>Medicon</b>	78.12.13
<b>SPECIFICATIONS</b>	<a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tendon Seizing Instruments</a>
<b>Aesculap</b>	BF307
<b>KLS Martin</b>	23-680-13
<b>Medicon</b>	78.12.13

## SPECIFICATIONS

[straight](#)