

TENDON INTERLACING FORCEPS

straight; double convertible

20.5cm (8.2")



SKU: 195-238-20

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Tendon Seizing Instruments |
| KLS Martin | 23-678-00 |
| Medicon | 78.14.20 |
| SPECIFICATIONS | straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Tendon Seizing Instruments |
| KLS Martin | 23-678-00 |
| Medicon | 78.14.20 |
| SPECIFICATIONS | straight |