

TENDON INTERLACING FORCEPS

gebogen; double convertible

20.5cm (8.2")



SKU: 195-239-20

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Seizing Instruments
KLS Martin	23-678-01
Medicon	78.14.21
SPECIFICATIONS	curved



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Seizing Instruments
KLS Martin	23-678-01
Medicon	78.14.21
SPECIFICATIONS	curved