

TENDON STRIPPER

Ø 5.0mm; flexible

23.0cm (9.2")

SKU: 195-265-05

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Stripper
KLS Martin	23-654-45
Medicon	78.20.55

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Stripper
KLS Martin	23-654-45
Medicon	78.20.55