

TENDON STRIPPER

Ø 5.0mm; flexible

23.0cm (9.2")

SKU: 195-265-05

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|---------------------------------|
| UNIT | pcs. |
| TYPE | Tendon Stripper |
| KLS Martin | 23-654-45 |
| Medicon | 78.20.55 |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|---------------------------------|
| UNIT | pcs. |
| TYPE | Tendon Stripper |
| KLS Martin | 23-654-45 |
| Medicon | 78.20.55 |