

## TENDON STRIPPER

Ø 4.0mm; flexible

23.0cm (9.2")

**SKU:** 195-265-04

**Categories:** [Open Surgery](#), [Orthopaedics](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tendon Stripper</a>
<b>KLS Martin</b>	23-654-40
<b>Medicon</b>	78.20.54

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tendon Stripper</a>
<b>KLS Martin</b>	23-654-40
<b>Medicon</b>	78.20.54