

## TOENNIS NEEDLE HOLDER

straight; serrated; without catch

18.0cm (7.2")

**SKU:** 050-362-18

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM510
<b>KLS Martin</b>	20-574-18
<b>AUTHORS</b>	<a href="#">TOENNIS</a>
<b>Medicon</b>	11.34.18
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM510
<b>KLS Martin</b>	20-574-18

## AUTHORS

[TOENNIS](#)

## Medicon

11.34.18

## SPECIFICATIONS

[serrated](#), [straight](#)