

## TONG TYPE ELECTRODE BIPOLAR; 3MM; INSERT

for bipolar coagulation forceps

628-101-00 and 628-105-00; 450mm

No CE mark! No labeling!



**SKU:** 628-140-54

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">HF-Instruments, bipolar</a>
<b>TYPE</b>	<a href="#">Electrodes</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">HF-Instruments, bipolar</a>
<b>TYPE</b>	<a href="#">Electrodes</a>