

TONG TYPE ELECTRODE BIPOLAR; 3MM; INSERT

for bipolar coagulation forceps

628-101-00 and 628-105-00; 330mm

No CE mark! No labeling!



SKU: 628-140-53

Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|-----------------------------------------|
| UNIT | pcs. |
| RANGE | HF-Instruments, bipolar |
| TYPE | Electrodes |

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|-----------------------------------------|
| UNIT | pcs. |
| RANGE | HF-Instruments, bipolar |
| TYPE | Electrodes |