

## TONSIL NEEDLE; CURVED

22 Gauge (0.7mm); 85mm

**SKU:** 084-201-85

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Guide Needles</a>
<b>KLS Martin</b>	18-641-20

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Guide Needles</a>
<b>KLS Martin</b>	18-641-20