

TOURNIQUET LIGATURE CATCHER

Ø2.5mm; 26.0cm


SKU: 141-028-25

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Catchers, Tourniquets
KLS Martin	24-550-25



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Catchers, Tourniquets
KLS Martin	24-550-25