

TOURNIQUET LIGATURE CATCHER

Ø3.0mm; 26.0cm

SKU: 141-028-30

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Catchers, Tourniquets
KLS Martin	24-550-30



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Catchers, Tourniquets
KLS Martin	24-550-30