

## TOURNIQUET LIGATURE CATCHER

Ø3.0mm; 26.0cm

**SKU:** 141-028-30

**Categories:** [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Ligature Catchers, Tourniquets](#)

**KLS Martin** 24-550-30



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Ligature Catchers, Tourniquets](#)

**KLS Martin** 24-550-30