

TOURNIQUET LIGATURE CATCHER

Ø2.0mm; 26.0cm



SKU: 141-028-20

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Ligature Catchers, Tourniquets](#)

KLS Martin 24-550-20



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Ligature Catchers, Tourniquets](#)

KLS Martin 24-550-20