

TOWEL FORCEPS FOR PAPER CLOTHS

curved; ball/socket

13.0cm (5.2")




SKU: 030-024-13

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Towel Forceps
Aesculap	BF464
KLS Martin	13-924-14
Medicon	17.55.23
SPECIFICATIONS	curved



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Towel Forceps
Aesculap	BF464
KLS Martin	13-924-14
Medicon	17.55.23

SPECIFICATIONS

[curved](#)