

## TOWEL FORCEPS FOR PAPER CLOTHS

curved; ball/socket

13.0cm (5.2")


**SKU:** 030-024-13

**Categories:** [Hemostats](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Towel Forceps</a>
<b>Aesculap</b>	BF464
<b>KLS Martin</b>	13-924-14
<b>Medicon</b>	17.55.23
<b>SPECIFICATIONS</b>	<a href="#">curved</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Towel Forceps</a>
<b>Aesculap</b>	BF464
<b>KLS Martin</b>	13-924-14
<b>Medicon</b>	17.55.23

## SPECIFICATIONS

[curved](#)