

## TOWEL FORCEPS FOR PAPER CLOTHS

curved; ball/socket

11.5cm (4.5")




**SKU:** 030-024-11

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Towel Forceps</a>
<b>Aesculap</b>	BF465
<b>KLS Martin</b>	13-924-11
<b>SPECIFICATIONS</b>	<a href="#">curved</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Towel Forceps</a>
<b>Aesculap</b>	BF465
<b>KLS Martin</b>	13-924-11
<b>SPECIFICATIONS</b>	<a href="#">curved</a>