

## TOWEL FORCEPS WITH TUBING CLAMP

curved; sharp

11.0cm (4.4")



**SKU:** 030-005-11

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Towel Forceps</a>
<b>Aesculap</b>	BF442
<b>KLS Martin</b>	13-913-11
<b>Medicon</b>	17.65.10
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">sharp</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Towel Forceps</a>
<b>Aesculap</b>	BF442
<b>KLS Martin</b>	13-913-11
<b>Medicon</b>	17.65.10

## SPECIFICATIONS

[curved](#), [sharp](#)