

## TOYNBEE EAR SPECULA

fig. 3

5.0mm



**SKU:** 111-035-03

**Categories:** [Open Surgery](#), [Otology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Ear Specula</a>
<b>Aesculap</b>	OF127
<b>KLS Martin</b>	36-118-12
<b>AUTHORS</b>	<a href="#">TOYNBEE</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Ear Specula</a>
<b>Aesculap</b>	OF127
<b>KLS Martin</b>	36-118-12
<b>AUTHORS</b>	<a href="#">TOYNBEE</a>