

TRIPOD POLYP RETRIEVER

teflon tube, with 3 sharp hooks

1.8mm, length 250cm; single use



SKU: PK-7GZEW1835250C-S

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Instruments for Double-Balloon-Enteroscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Instruments for Double-Balloon-Enteroscopy](#)