

TRIPOD POLYP RETRIEVER;

teflon tube, with 8 sharp hooks

2.6mm, length: 230cm

with solid fix click handle



SKU: PK-7GZRG26252300KW-Z

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Foreign Body Forceps](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Foreign Body Forceps](#)