

## TRIPOD POLYP RETRIEVER;

teflon tube, with 8 sharp hooks

2.6mm, length: 230cm

with solid fix click handle



**SKU:** PK-7GZRG2625230OKW-Z

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Foreign Body Forceps](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Foreign Body Forceps](#)