

## TRIPOD POLYP RETRIEVER (SINGLE USE)

with 4 round hooks

tubus 2.5 x 230cm



**SKU:** PK-7GREW2540230Q-S

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)