

## TRIPOD POLYP RETRIEVER WITHOUT HANDLE;

teflon tube, with 4 sharp hooks

2.5mm, length: 230cm



**SKU:** PK-7GZPK2540230Q

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Foreign Body Forceps](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Foreign Body Forceps](#)