

TRIPOD POLYP RETRIEVER WITHOUT HANDLE;

teflon tube, with 4 sharp hooks

2.5mm, length: 230cm



SKU: PK-7GZPK2540230Q

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Foreign Body Forceps](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Foreign Body Forceps](#)