

## TRIPOD POLYP RETRIEVER WITHOUT HANDLE;

teflon tube, with 3 round hooks

1.8mm, length: 120cm



**SKU:** PK-7GRPK1840120C

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Foreign Body Forceps](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Foreign Body Forceps](#)