

## TUBING CLAMP

straight; serrated

18.0cm (7.2")



**SKU:** 032-002-18

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)


**TYPE** [Tubing Clamps](#)

**Aesculap** AN912R

**KLS Martin** 13-944-18

**Medicon** 17.72.78

**SPECIFICATIONS** [serrated](#), [straight](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Tubing Clamps](#)

**Aesculap** AN912R

**KLS Martin** 13-944-18

**Medicon** 17.72.78

## SPECIFICATIONS

[serrated](#), [straight](#)