

TUBING CLAMP

straight; serrated

20.0cm (8")




SKU: 032-002-20

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Tubing Clamps |
| Aesculap | AN913 |
| KLS Martin | 13-944-20 |
| Medicon | 17.72.80 |
| SPECIFICATIONS | serrated , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-------------------------------|
| UNIT | pcs. |
| TYPE | Tubing Clamps |
| Aesculap | AN913 |
| KLS Martin | 13-944-20 |
| Medicon | 17.72.80 |

SPECIFICATIONS

[serrated](#), [straight](#)