

TUBING CLAMP

straight; serrated

20.0cm (8")




SKU: 032-002-20

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
Aesculap	AN913
KLS Martin	13-944-20
Medicon	17.72.80
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
Aesculap	AN913
KLS Martin	13-944-20
Medicon	17.72.80

SPECIFICATIONS

[serrated](#), [straight](#)