

TUBING CLAMP

straight; serrated; with Clip

15.5cm (6.2")



SKU: 032-004-15

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Tubing Clamps](#)

KLS Martin 13-946-15

Medicon 17.72.76

SPECIFICATIONS [serrated](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Tubing Clamps](#)

KLS Martin 13-946-15

Medicon 17.72.76

SPECIFICATIONS [serrated](#), [straight](#)