

## TUBING CLAMP

straight; serrated; with Clip

20.0cm (8")



**SKU:** 032-004-20

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN


**UNIT** [pcs.](#)

**TYPE** [Tubing Clamps](#)

**KLS Martin** 13-946-20

**Medicon** 17.72.81

**SPECIFICATIONS** [serrated](#), [straight](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Tubing Clamps](#)

**KLS Martin** 13-946-20

**Medicon** 17.72.81

**SPECIFICATIONS** [serrated](#), [straight](#)