

TUBING CLAMP

straight; serrated; with Clip

20.0cm (8")



SKU: 032-004-20

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Tubing Clamps](#)

KLS Martin 13-946-20

Medicon 17.72.81

SPECIFICATIONS [serrated](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Tubing Clamps](#)

KLS Martin 13-946-20

Medicon 17.72.81

SPECIFICATIONS [serrated](#), [straight](#)