

TUBING CLAMP

straight; serrated

15.5cm (6.2")



SKU: 032-002-15

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
Aesculap	AN910
KLS Martin	13-944-15
Medicon	17.72.75
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
Aesculap	AN910
KLS Martin	13-944-15
Medicon	17.72.75

SPECIFICATIONS

[serrated](#), [straight](#)