

## TUOHY BORST ADAPTER

with LuerLock connection (stainless steel)

incl. 3x Silicone seals



**SKU:** 630-047-30

**Categories:** [Endoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Instruments for Laser Therapy</a>
<b>TYPE</b>	<a href="#">Working Elements</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Instruments for Laser Therapy</a>
<b>TYPE</b>	<a href="#">Working Elements</a>