

## TWIST DRILL

smooth shaft

Ø 6.0mm

**SKU:** 190-288-60

**Categories:** [Open Surgery](#), [Orthopaedics](#)



## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hand Drills</a>
<b>KLS Martin</b>	22-199-60
<b>Medicon</b>	76.69.60

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hand Drills</a>
<b>KLS Martin</b>	22-199-60
<b>Medicon</b>	76.69.60