

TWIST DRILL

smooth shaft

Ø 2.5mm

SKU: 190-288-25

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-----------------------------|
| UNIT | pcs. |
| TYPE | Hand Drills |
| KLS Martin | 22-199-25 |
| Medicon | 76.69.25 |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-----------------------------|
| UNIT | pcs. |
| TYPE | Hand Drills |
| KLS Martin | 22-199-25 |
| Medicon | 76.69.25 |