

TWIST DRILL

smooth shaft

Ø 2.8mm

SKU: 190-288-28

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-----------------------------|
| UNIT | pcs. |
| TYPE | Hand Drills |
| KLS Martin | 22-199-28 |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-----------------------------|
| UNIT | pcs. |
| TYPE | Hand Drills |
| KLS Martin | 22-199-28 |