

TWIST DRILL FOR HAND DRILL

Ø 2.0mm

SKU: 190-261-20

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
Aesculap	FR040
KLS Martin	22-145-12

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
Aesculap	FR040
KLS Martin	22-145-12