

## TWIST DRILL FOR HAND DRILL

Ø 2.0mm

**SKU:** 190-261-20

**Categories:** [Open Surgery](#), [Orthopaedics](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hand Drills</a>
<b>Aesculap</b>	FR040
<b>KLS Martin</b>	22-145-12

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hand Drills</a>
<b>Aesculap</b>	FR040
<b>KLS Martin</b>	22-145-12