

TWIST DRILL FOR HAND DRILL

Ø 3.0mm

SKU: 190-261-30

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
Aesculap	FR042
KLS Martin	22-145-02
Medicon	76.50.12

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
Aesculap	FR042
KLS Martin	22-145-02
Medicon	76.50.12