

## TWIST DRILL FOR HAND DRILL

Ø 3.5mm

**SKU:** 190-261-35

**Categories:** [Open Surgery](#), [Orthopaedics](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hand Drills</a>
<b>Aesculap</b>	FR043
<b>KLS Martin</b>	22-145-03
<b>Medicon</b>	76.50.13

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hand Drills</a>
<b>Aesculap</b>	FR043
<b>KLS Martin</b>	22-145-03
<b>Medicon</b>	76.50.13