

TWIST DRILL FOR HAND DRILL

Ø 3.5mm

SKU: 190-261-35

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
Aesculap	FR043
KLS Martin	22-145-03
Medicon	76.50.13

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
Aesculap	FR043
KLS Martin	22-145-03
Medicon	76.50.13