

## ULRICH SPONGE FORCEPS

straight

27.0cm (10.8")


**SKU:** 030-200-27

**Categories:** [Hemostats](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Sponge and Dressing Forceps</a>
<b>Aesculap</b>	BF062
<b>KLS Martin</b>	14-100-27
<b>AUTHORS</b>	<a href="#">ULRICH</a>
<b>Medicon</b>	19.40.26
<b>SPECIFICATIONS</b>	<a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Sponge and Dressing Forceps</a>
<b>Aesculap</b>	BF062
<b>KLS Martin</b>	14-100-27

## AUTHORS

[ULRICH](#)

## Medicon

19.40.26

## SPECIFICATIONS

[straight](#)