

URETER SLITTING SHEAR; Ø 1.7MM; FL 600M

5 Fr.; upwards cutting; semi-rigid

single action; with LuerLock;with handle



SKU: 583-490-16LL

Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Uretero-Renoscopy
TYPE	Semi-rigid Instruments



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Uretero-Renoscopy
TYPE	Semi-rigid Instruments