

## URETER SLITTING SHEAR; Ø 1.7MM; FL 600M

5 Fr.; upwards cutting; semi-rigid

single action; with LuerLock;with handle



**SKU:** 583-490-16LL

**Categories:** [Endoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Uretero-Renoscapy](#)

**TYPE** [Semi-rigid Instruments](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Uretero-Renoscapy](#)

**TYPE** [Semi-rigid Instruments](#)